



medical systems

From Azenta Life Sciences

B Connected

24/7 real-time monitoring solution for the safety of samples



- 24/7 real-time monitoring of internal temperature, door status, alarms, and warnings
- Monitor temperatures ranging from -99°C to +99°C¹
- Web-based portal.² Compatible with web-based and mobile browsers
- Instant alerts via SMS and email
- Easy data access and export in the form of .pdf, .csv, .xml and .xlsx reports
- Raw data access via REST API
- Unlimited data storage and history³
- Audit trail functionality
- Also compatible with non-B Medical Systems devices⁴
- Available in English, French, Spanish, German, Russian and Czech

Developed according to
GAMP® 5

21 CFR Part 11 Subpart
A and B Compliance

¹. Based on the sensors
². Installed on a server provided by the customer
³. Depends on the data storage capacity of the customer's server
⁴. Via temperature sensors connected to a Data Communication Unit. Conditions apply

SAVING LIVES THROUGH RELIABLE AND INNOVATIVE TECHNOLOGY

B Medical Systems is a global manufacturer of medical grade refrigeration and transport solutions with over **40 years of experience.**



B Medical Systems S.à r.l

17, op der Hei, L-9809, Hosingen, Luxembourg



www.bmedicalsistemas.com